

# ISOLUTE® QuEChERS Method

## EN 15662 Highly Pigmented Fruit and Vegetables

### This product contains:

PSA	150 mg
MgSO <sub>4</sub>	900 mg
GCB	45 mg

### Part numbers required:

Q0020-15V	10 g QuEChERS EN 15 mL Extraction Tube
Q0090-15V	EN Highly Pigmented Fruit and Vegetables Clean Up Tube

### To be supplied by the User:

- Internal standard
- Acetonitrile
- Centrifuge
- Vortex
- 50 mL centrifuge tubes (part number Q0000-50V) or equivalent



### ISOLUTE QuEChERS Method

## EN 15662 Highly Pigmented Fruits and Vegetables

Part Numbers: Q0090-15V and Q0020-15V

1	Macerate food sample in a suitable container to a smooth paste.
2	Weigh 10 g of macerated sample into an empty 50 mL centrifuge tube or other suitable container
3	Add internal standard
4	Add 10 mL of Acetonitrile
5	Pour the contents of the 15 mL tube (Q0020-15V) onto the sample
6	Shake by hand for 1 minute
7	Centrifuge 4000 rpm for 5 minutes
8	Transfer 6 mL of supernatant to the 15 mL dSPE tube (Q0090-15V)
9	Vortex 1 minute
10	Centrifuge 4000 rpm for 5 minutes
11	Transfer aliquot of supernatant to GC vial for analysis